

# VETERINARY INSTRUMENTATION

## SUMMER SPECIAL OFFERS

### CrCL Suture System



Other kits available.

The lateral suture remains the technique of choice for most veterinary surgeons worldwide.

#### Contains:

- Crimping Forceps
- Heavy Duty Needleholders
- 2 of each size (50, 80, 100lb) sterile nylon leader line pack with swaged on needle and crimp.
- Procedural DVD

Code: TA091151

Components Price	£263.00
Kit Price	£170.00
<b>OFFER PRICE</b>	<b>£125.00</b>

### Sterile CrCL Multipacks



Code: TA091159

Each Multipack contains the following:

- 5 x 50lb nylon line x 500mm on swaged-on v.small fabella needle +10mm crimp sterile
- 5 x 80lb nylon line x 800mm on swaged-on small fabella needle +12mm crimp sterile
- 5 x 100lb nylon line x 800mm on swaged-on medium fabella needle +12mm crimp sterile

Each pack contains everything for a single procedure.

Components Price	£232.50
Kit Price	£150.00
<b>OFFER PRICE</b>	<b>£125.00</b>

### LigaFiba Lateral Button Suture Kit



This new lateral suture technique has been developed to minimise the tendency for the proximal strand of the suture to drift up the straight patella ligament. Step by step instruction sheet included.

#### Contains:

- LigaFiba Crimping Forceps
- 2x each size LigaFiba /needle /crimp /button sterile pack (150lb, 250lb, 500lb).
- Step by step instruction sheet

Components Price	£320.00
Kit Price	£250.00

**OFFER PRICE £175.00**

Code: TALFLBSKIT

### LigaFiba Isotoggle Starter Kit



The LigaFiba Isotoggle suture is placed via bone tunnels across the stifle at isometric points. The use of isometric points maximises the range of movement and protects the suture.

#### The Starter Kit contains:

- 2 x 250lb Isotoggle sutures
- 2x 500lb Isotoggle sutures
- 1x 2.5mm antiskid drill
- 1x 3.5mm antiskid drill

Components Price	£347.00
Kit Price	£275.00

**OFFER PRICE £225.00**

Code: TALFISOKIT

**LIFETIME WARRANTY AGAINST MANUFACTURING DEFECT ON ALL INSTRUMENTS. FREE DELIVERY (UK)**

**VALID UNTIL 30/09/11**

**Tel: 0845 130 9596 Fax: 0845 130 8687 E-mail: info@vetinst.com Website: www.vetinst.com**

Please note offers cannot be used in conjunction with any other offers or discounts. Offers only available if ordered directly rather than through wholesalers. All prices exclude VAT

# SUMMER SPECIAL OFFERS

VALID UNTIL 30/09/11 [www.vetinst.com](http://www.vetinst.com) Tel: 0845 130 9596

## Suture Anchor Set



Suture screws are a versatile, economic suture anchor. Used together with LigaFiba, suture screws offer a solution to joint stabilisation. **Set contains:**

- 2 x 2.0mm diameter 6 & 10mm long
- 2 x 2.7mm diameter 8 & 14mm long
- 2 x 3.5mm diameter 12 & 20mm long
- 2 x 4.00mm diameter 16 & 24mm long
- Universal Suture Screw Introducer
- LigaFiba sterile suture packs x 1 of each size (150lb, 250lb, 500lb)



**Code: TASSITKIT**

Components Price £510.00  
Kit Price £395.00

**OFFER PRICE £335.00**

## Meniscus Surgery Kit



**Contains:**

- Standard Stifle Distractor
- Meniscus Probe 1mm
- Beaver handle 130mm
- 8mm Hohman retractor
- 1 Halstead forceps 1/2 teeth, straight
- Stainless box
- Silicone insert & drape
- 25 x SM65 blades (mini no.11)

A complete Meniscus Surgery Kit containing purpose designed instrumentation, all conveniently contained in a stainless steel case with protective silicone lining.

**Code: TA001116**

Components Price £390.75  
Kit Price £295.00

**OFFER PRICE £250.00**

## External Stifle Distractor



Osteoarthritic stifles can be very 'tight' making examination and treatment of meniscal lesions difficult. Distraction of the articular surfaces using intra-articular instruments clutter the operative field.

Pins are placed above and below the stifle approximately along the line of the medial collateral ligament. The distractor uses 2.5mm pins (Max).

**Code: TAAR-17-01**

Component Price £145.00

**OFFER PRICE £105.00**

## Bone Plate / Drill Guide Set



Temporary fixation of a bone plate at a fracture is an essential part of bone plating. Unfortunately there is an equal desire for the placement of bone reduction forceps and drill guides. The new range of VI bone holding / drill guide forceps attempts to address two of these requirements.

**Contains:**

- 1.5/2.0 Plate Holding / Drill Guide Forceps
- 2.4/2.7 Plate Holding / Drill Guide Forceps
- 3.5 Plate Holding / Drill Guide Forceps

**Code: TA0012SET**

Components Price £315.00  
Kit Price £250.00

**OFFER PRICE £195.00**

## Cruciate Starter Kits

**TTO**



**TPLO**



**TTA**




CODE	DESCRIPTION	COMPONENT PRICE	KIT PRICE	SPECIAL OFFER PRICE
<b>TATTOBKIT</b>	TTO Basic Starter Kit	£1,040.00	£875.00	<b>£775.00</b>
<b>TATTOBKITPOW</b>	TTO Basic Starter Kit plus Multisaw Kit	£1,515.00	£1,350.00	<b>£1,175.00</b>
<b>TATPLODKIT</b>	TPLO Starter Kit with Delta Plates	£4,509.50	£3,750.00	<b>£3,500.00</b>
<b>TATPLOSKIT</b>	TPLO Starter Kit with Slocum Plates	£4,417.50	£3,750.00	<b>£3,500.00</b>
<b>TATTAPREM</b>	TTA Starter Kit in Premium Box contains all instruments/implants	£4,091.75	£3,475.00	<b>£3,175.00</b>

Visit [www.vetinst.com](http://www.vetinst.com) for complete kit lists

# SUMMER SPECIAL OFFERS

VALID UNTIL 30/09/11 [www.vetinst.com](http://www.vetinst.com) Tel: 0845 130 9596

## Locking Plate and Screw Set



The locking plate and screw set contains a wide selection of locking plates and screws as well as all the essential instruments required. Available in both 2.7mm and 3.5mm. **3.5mm set contains:**

- 1 of each plate (4 - 14 hole)
- 3 of each screw 12mm to 20mm
- 4 of each screw 22mm to 26mm
- 5 of each screw 28mm to 38mm
- 4 of each screw 40mm to 48mm
- 3 of each screw 50mm to 58mm
- 5 locking plugs
- 1x 2.5mm Hex screwdriver
- 2.7mm Hard Drill
- Drill guide
- 1x Screw holding forceps
- Premium Box / Insert

### 3.5 Locking Plate and Screw Set

Code: TALS35KITALL

Components Price £2,438.00

Kit Price £2,150.00

**OFFER PRICE £2,000.00**

### 2.7 Locking Plate and Screw Set

Code: TALS27KITALL

Components Price £1,941.25

Kit Price £1,725.00

**OFFER PRICE £1,500.00**

## Cortical Plating Kits



All kits are available as either cortical or cortical self tapping. The description cortical means that the pilot hole must be pre-tapped prior to screw placement. Cortical self tapping means that the screw itself has a self tapping tip and cuts its own thread. Self tappers work very well in all but the hardest cortical bone. It is recommended that surgeons using self tappers keep a tap 'just in case.'

Each kit contains all you need including screws and implants.

CODE	DESCRIPTION	COMPONENT PRICE	KIT PRICE	SPECIAL OFFER PRICE
TACORTKIT1.5/2.0	1.5/2.0mm Cortical Kit in Stainless Box	£1,728.00	£1,625.00	<b>£1,425.00</b>
TACORTKIT2.0	2.0mm Cortical Kit in Stainless Box	£1,310.00	£1,200.00	<b>£1,050.00</b>
TACORTKIT2.7/3.5	2.7/3.5mm Cortical Kit in Stainless Box	£2,261.00	£2,100.00	<b>£1,925.00</b>
TACORTKIT4.5	4.5mm Cortical Kit in Stainless Box	£1,352.00	£1,250.00	<b>£1,075.00</b>
TACORTKIT4.5/5.5	4.5/5.5mm Cortical Kit in Stainless Box	£2,602.00	£2,425.00	<b>£2,200.00</b>
TACORTSTKIT1.5/2.0	1.5/2.0mm Cortical Self Tapping Kit in Stainless Box	£1,635.00	£1,425.00	<b>£1,275.00</b>
TACORTSTKIT2.0	2.0mm Cortical Self Tapping Kit in Stainless Box	£1,215.00	£1,100.00	<b>£975.00</b>
TACORTSTKIT2.4	2.4mm Cortical Self Tapping Kit in Stainless Box	£1,344.00	£1,250.00	<b>£1,075.00</b>
TACORTSTKIT2735	2.7/3.5mm Cortical Self Tapping Kit in Stainless Box	£2,047.00	£1,925.00	<b>£1,750.00</b>
TACORTSTKIT4.5	4.5mm Cortical Self Tapping Kit in Stainless Box	£1,255.00	£1,100.00	<b>£975.00</b>

Visit [www.vetinst.com](http://www.vetinst.com) for complete kit lists

## MiniDriver Sets



### The MiniDriver is Dead! Long Live the Minidriver!

The entire MiniDriver range is now available and when purchased as a set the savings are considerable.

The entire range is remanufactured and comes with a 12 month warranty.

#### Set Contains:

Handpiece, 4mm chuck, wire driver, AO quick release chuck, Sagittal saw Mark II, Minidriver Hose.

Code: TAVIMDB2

Components Price £5035.00

Kit Price £3750.00

**OFFER PRICE £3395.00**

## Blue Line Battery Surgical Power



The Blue Line drill and saw have proved very popular since their launch in 2010.

The drill is fully autoclavable and is supplied with a quick fit 1/4" (6.3mm) stainless chuck and an AO quick release chuck. Fully cannulated.

The saw accepts a range of blades from 6mm wide. Saw speed is nicely controlled via the trigger.

Supplied with 12 month warranty.

#### Set Contains:

Drill, Saw, Jacobs chuck, AO quick release chuck, 2x battery chargers, 3x batteries 2x battery transfer sleeves.

Code: TABLDSSET

Components Price £3415.00


Kit Price £2750.00

**OFFER PRICE £2450.00**

# SUMMER SPECIAL OFFERS

VALID UNTIL 30/09/11 [www.vetinst.com](http://www.vetinst.com) Tel: 0845 130 9596

## Whiting Bone Plate Retractor Set



The Whiting Bone Plate Retractor is designed to expose the end hole(s) of the plate for screw placement. It is often the case that a plate is pushed under soft tissues to avoid unnecessarily large exposures only to find that it is very difficult to access the last screw hole. This retractor lifts the soft tissue to allow access.

**Code: TABPRSET**

### Set Contains:

Retractor for 3.5 / 3.5 Broad  
Retractor for 2.4/2.7

Components Price £125.00

**OFFER PRICE £105.00**

## Orthopaedic Wire Kit



### Contains:

- Wire twister
- Set of 3 single ended wire passers
- 10m roll of 1.0mm orthopaedic wire

**Code: TAOWKIT**

Cerclage wire can be a very useful means of fragment fixation.


This kit contains our most popular cerclage instruments and a roll of wire.

Components Price £178.00

Kit Price £135.00

**OFFER PRICE £105.00**

## Strong Wire Passer Set



This three instrument set is double ended giving internal diameters of 12.5 & 15mm on the smallest, 20 and 25mm (medium) and 30 and 40mm (Large).

**Code: TACWPSET**

Components Price £178.00

Kit Price £125.00

**OFFER PRICE £105.00**

## Radiographic Reference Ball



Taking radiographs where the area of interest is not close to the plate creates some challenges in scaling. This is true of both digital and analogue systems. The Reference Ball is 25mm in diameter and is attached to a stable base by a flexible neck which allows the ball to be placed at the same level above the plate as the area of interest.

Using a spherical object as a reference eliminates errors created by ruler type reference which might not be quite parallel with the plate.

**Code: TAXRAY25MM**

Components Price £60.00

**OFFER PRICE £50.00**

## Safety Pin Cutters - New



The new telescopic handles make it easier to use, easier to clean and easier to autoclave.

Same cutting action as original.

**Code: TA001229**

Components Price £225.00

**OFFER PRICE £175.00**

## Modular Osteotome



The osteotome is supplied as a comfortable handle into which fits a large range of blades. The blades are very thin 1.1mm making them suitable for a variety of procedures including block sulcoplasty.

Supplied as a handle with the following blades 4, 6, 8, 10, 12, 15 and 20mm.

A diamond-coated rasp and 2 & 3mm blades also available.

**Code: TA001380**

Components Price £170.00

**OFFER PRICE £135.00**

# SUMMER SPECIAL OFFERS

VALID UNTIL 30/09/11 [www.vetinst.com](http://www.vetinst.com) Tel: 0845 130 9596

## Aural Haematoma System



Surgical management of aural haematomas is a common procedure in veterinary practice. Once drained the haematoma must be prevented from filling up again. This new system developed by Nick Hart makes the procedure straight forward.

The system comprises a pair of clamp forceps and 24 implants.

The forceps apply a fenestrated plastic implant under controlled pressure which is then stitched into place.




Code: TA014600

Components Price £163.00  
Kit Price £115.00

**OFFER PRICE £105.00**

## Ring Retractor Kit



Some abdominal approaches require a more complex retraction system than a simple Balfour. The VI Ring Retractor allows the surgeon to place hooked retractor strips around the periphery of a surgical incision which may be tensioned individually and retained around the edge.

Code: TA701700

Contains:

- 3 Standard rings ( 200 x 130mm, 160 x 85mm, 110 x 40mm internal)
- Set (10) sterile hooks

Components Price £167.50  
Kit Price £135.00

**OFFER PRICE £105.00**

## AHM Balfour Retractor



Alasdair Hotson Moore has re-visited the original Balfour design and has made some significant changes which improve the instrument both in use and sterilisation.

Code: TAAHM

Components Price £155.00

**OFFER PRICE £125.00**

## Soft Edged Surgical Spoon Set



Removing calculi from the bladder or impacted faeces from the rectum requires great care to avoid damage. 'Domestic' spoons are typically too sharp and too short. Available as set of 2 small (20x25x250mm) and large (67x34x270mm).

Code: TA014702

Components Price £39.00  
Kit Price £30.00

**OFFER PRICE £25.00**

## Comfort Handle Elevators



Designed and developed by David Strong. The comfort pattern elevators takes a radical new look at elevator handles. The pattern fits into the palm of the hand allowing the fore finger to be positioned along the shaft for added control. The kit contains 1, 2, 3, 4mm widths in a stainless case.

Code: TA012600

Components Price £210.00  
Kit Price £150.00

**OFFER PRICE £105.00**

## Periodontal Hand Instrument Set



Seven periodontal instruments in a stainless autoclave case.

Instruments include mirror, shepherd hook explorer/measuring probe, straight curette scaler, angled curette scaler, angled jaquette scaler plus two fine probes.

Code: TA012210

Components Price £135.00

**OFFER PRICE £105.00**

## Spay Kits



These comprehensive kits contain all necessary instrumentation.

### Feline Spay Kit

- Stainless Autoclave Case
- Metz Scissor Str 14.5cm
- Sturdy type spay forceps
- Adson rat tooth forceps
- Halstead Mosquito Forceps x 2
- Fenestrated drape
- Cross action towel clips x 2
- Olsen Hegars 15.5cm
- Spencer Wells Curved x 2
- No.3 Scalpel handle

Code: TASPAYFEL

Components Price	£294.00
Kit Price	£255.00

**OFFER PRICE £185.00**

### Canine Spay Kit

- Stainless Autoclave Case
- Curved Hysterectomy Clamps x 2
- Mayo Scissors 14cm Straight
- Rat Toothed Treves 14cm
- Spencer Wells 14cm Curved x 4
- Allis Tissue Forceps x 2
- Fenestrated Drape
- Olsen Hegars 16.5cm
- Cross Action Towel Clips x 4
- No.3 Scalpel handle

Code: TASPAYCAN

Components Price	£445.00
Kit Price	£375.00

**OFFER PRICE £285.00**

**BUY BOTH THE FELINE AND CANINE SPAY KITS AND SAVE AN ADDITIONAL 10%**  
Code: TASPAYBOTH

## Ophthalmic Kit



### Contains:

- Roberts Tying Forceps
- Castroviejo corneal scissors straight
- Halstead mosquito forceps curved and straight
- Barraquer wire eyelid speculum
- Beaver type handle
- Castroviejo needleholders with catch straight
- Steven's tenotomy scissors straight
- Colibri tissue forceps
- Graefe fixation forceps 3.5mm jaws
- Stainless case
- Silicone protective insert

The kit contains 10 instruments. All are contained in a stainless steel case with a protective silicone lining which minimises damage to the delicate instruments. As recommended by David Williams MA VetMB PhD CertVOphthal, for a range of ophthalmic procedures.

All instruments have a lifetime warranty against manufacturing defect.

Components Price	£698.00
Kit Price	£455.00

**OFFER PRICE £365.00**

Code: TA014130

## Rodent Dental Kit



With rodents increasing in popularity, this kit contains all the instruments necessary for rodent dentistry. All instruments have a lifetime warranty against manufacturing defect.

### Contains:


- Small fine diamond molar rasp
- Rodent spatula
- Cheek dilators small, large & long reach
- Crossley incisor & molar luxators
- Molar extraction forceps & cutters
- Mouth gag
- Instrument box

Code: TAROD2

Components Price	£335.00
Kit Price	£250.00

**OFFER PRICE £185.00**

## Economy Table Top Rodent Gag & Light



Allows rodent mouth to be held open and illuminated in an appropriate position without assistance. White teflon platform. Table fully autoclavable.

**BUY BOTH THE RODENT DENTISTRY KIT AND TABLE TOP GAG & LIGHT AND SAVE AN ADDITIONAL 10%**

Code: TAROBTOTH

Code: TA012209

Components Price	£305.00
Kit Price	£275.00

**OFFER PRICE £235.00**

## VETERINARY INSTRUMENTATION

Broadfield Road,  
Sheffield

S8 0XL U.K.

Tel: 0845 130 9596

Fax: 0845 130 8687

Overseas Customers use

Tel: +44 114 258 8530

Fax: +44 114 255 4061

E-mail: [info@vetinst.com](mailto:info@vetinst.com)

Website: [www.vetinst.com](http://www.vetinst.com)